

MAXcyli

MEDICAL

ALUMINUM OXYGEN CYLINDERS

Aluminum cylinders have become the cylinder of choice in the medical industry due to their clean finish, lightweight and easy to handle characteristics. ProRack Gas Control Products offers many sizes with various valves depending on the customers preference. When you are in the market for cylinders, CALL ProRack TODAY!

- Large Inventory
- Quick Delivery
- Great Prices
- Fast Shipping Nationwide!

CYLINDER NOTES

- All cylinders sized from M2 to E are designed with a .750-16UNF2A inlet thread.
- All cylinders sized from M60 to MM are designed with a 1.125-12UNF2A inlet thread.
- All cylinders are supplied standard with a protective clear coat with a green shoulder.
- All cylinders are DOT/TC compliant and stamped accordingly for the intended service.

CYLINDER OPTIONS

- Custom powder coat painting is available.
- Stamping is available.
- Label under the clear coat is available.
- Customizing available.

VALVE OPTIONS

- Copper siphon tubes are available at additional fee.
- Standard brass valves from Sherwood and Cavagna are available. CDI stocks a wide selection of common valves such as CGA320, 326, 346, 350, 540, 580, 590 & 660.
- RPV valves, brass diaphragm valves and stainless steel diaphragm valves are available.





ALUMINUM OXYGEN CYLINDER ORDERING INFORMATION *

CGA870 Post Valve Options							CGA540 Valve Options	
Cylinder	Wrenchtight	Toggle	Handtight	Gauge	Gauge & Toggle	Gauge & Handtight	Low Profile w/gauge	Standard
M2	K870M2	K870M2T	K870M2H	K870M2G	K870M2GT	K870M2GH	Y540M2G	T540M2
M4	K870M4	K870M4T	K870M4H	K870M4G	K870M4GT	K870M4GH	Y540M4G	T540M4
M6	K870M6	K870M6T	K870M6H	K870M6G	K870M6GT	K870M6GH	Y540M6G	T540M6
ML6	K870ML6	K870ML6T	K870ML6H	K870ML6G	K870ML6GT	K870ML6GH	Y540ML6G	T540ML6
M7	K870M7	K870M7T	K870M7H	K870M7G	K870M7GT	K870M7GH	Y540M7G	T540M7
M9	K870M9	K870M9T	K870M9H	K870M9G	K870M9GT	K870M9GH	Y540M9G	T540M9
D	K870D	K870DT	K870DH	K870DG	K870DGT	K870DGH	Y540DG	T540D
Jumbo D	K870DJ	K870DJT	K870DJH	K870DJG	K870DJGT	K870DJGH	Y540DJ	T540DJ
E	K870E	K870ET	K870EH	K870EG	K870EGT	K870EGH	Y540EG	T540E
M60	-	-	-	-	-	-	-	M60-540
MM	-	-	-	-	-	-	-	MM540

* All part numbers shown are cylinders with oxygen valves installed.

ALUMINUM OXYGEN CYLINDER SPECIFICATIONS WITHOUT VALVES

Cylinder	Skid Qty	Oxygen Capacity		Service PSI		Cylinder Length		Cylinder O.D.		Cylinder Weight	
		CU.FT.	LITER	PSI	BAR	IN.	MM	IN.	MM	LBS	KG
M2	20 per box	1.5	42	2216	153	5.6	143	2.5	63.5	.74	.34
M4	144	4	113	2216	153	8.4	213	3.2	81	1.6	7.4
M6	144	6	170	2216	153	11.59	294	3.2	81	2.2	1
ML6	81	6	170	2015	139	7.68	195	4.38	111	2.9	1.29
M7	81	7	198	2015	139	9.18	233	4.38	111	2.9	1.29
M9	81	9	255	2015	139	11.88	276	4.38	111	3.7	1.69
D	81	15	425	2015	139	16.51	419	4.38	111	5.3	2.41
Jumbo D	72	22.6	623	2015	139	16.3	414	5.25	133	8.1	3.68
E	81	24	680	2015	153	25.63	651	4.38	111	7.9	3.58
M60	25	60	1699	2216	153	23	584	7.25	184	21.7	9.86
MM	25	122	3455	2216	153	35.75	908	8	203	38.6	17.55

* All capacity specifications are based on oxygen service.

Santa Ana, CA
Sacramento, CA
Baltimore, MD

1300 East Wakeham, Unit A
Santa Ana, CA 92703

Gas and Cryogenics **800.752.4579**
www.prorackgasproducts.com